

BrailleNote Touch Plus 2.6 – Selecting, Deleting, and Formatting Text

Selecting Text			
	SELECTING TEXTNavigate your cursor to the BEGINNING of the portion of text you would like to select. (see Tutorial 2.5 for instructions on Navigation within a document.)Press ENTER + S to activate selection mode.Navigate your cursor to the END of the text you would like to select. (for example, Spacebar + 6 moves you to the next word.)Press ENTER + S again to finish your selection.You will now see DOTS 7 and 8 on the refreshable braille display "flashing" below your selected text.	HumaRuare Balakoker fordh Image: Comparison of the second seco	
	SELECT ALL To select all text in a document, press ENTER + DOT 123456.		

Deleting Text DELETE SELECTED TEXT Select a portion of text, then press BACKSPACE.

SET►BC

DELETE LETTER To delete the letter to the left of your cursor, press BACKSPACE.	
DELETE CURRENT CHARACTER To delete the current character, press BACKSPACE + 36.	
DELETE PREVIOUS WORD To delete the word to the left of your cursor, press BACKSPACE + 2.	
DELETE CURRENT WORD To delete the word your cursor is on, press BACKSPACE + 25.	
DELETE REMAINDER OF LINE To delete text from your cursor to the end of the line, press BACKSPACE + 14 .	
DELETE REMAINDER OF DOCUMENT To delete text, from your cursor to the end of the document, press BACKSPACE + DOT 456. Note: Before deleting the remainder of the document, you will get a message asking if you are sure you want to delete to the end of the document. Navigate to "yes" to complete the action.	

Formatting Text


FORMATTING TEXT

When a portion of text is highlighted, you can change the format of that text by pressing one of the following commands: press ENTER + B for bold text

press ENTER + i for italicized text

press ENTER + U for underlined text

To access all the options for formatting, you can access the Context Menu>Format Functions by pressing **SPACE + M**, then navigate to "Format Functions." (Or long press the square Recent Apps button to go to the same place.)



Tip...

The Humanware Manual for the BNT+ can be found here, including a full list of shortcuts. http://support.humanware.com/Site/Files/a/06cc20404d351cbaa2520fcd0d5adc3/f117ef 24d066dbd8a70fcaf08af3d539/EN-UG%20BNTouchPlus.pdf