

Communication Device/Hardware	Operating System	Sizes Available	Communication Software available	Access options	Environmental Controls capability	Keyguards/Touchguide available
Accent (prentke romich)	Windows 10	8 inch/800 10 inch/1000 14 inch/1400	Nuvoice: which includes a large number of communication page sets.	Touch/mouse/joystick 1 and 2 switch access Head mouse-Nupoint (add on) Eye gaze (Nueye/add on) Auditory prompts, Auditory scanning	yes	Yes Key guards and touch guides available
Tobii Dynavox I series	Windows 10	12 inch/ I-12 15 inch/ I-15	Communicator 5 Grid 3 ?? Compass or Compass w/ PODD Snap + Core First	Touch/mouse/joystick, 1 and 2 switch access, head mouse (add on), eye gaze (add-on)(only windows control classic software available with eye gaze)	yes	no
NovaChat (Saltillo)	Android	5, 8, 10, 12	NovaChat app LAMP App	Touch, 1 and 2 switch access, auditory scanning	no	yes
iPad	iOS	-iPad mini/ 7.9 inch diag. -iPad air /9.7 inch -iPad pro/ 12.9 inch	LAMP, TouchChat, Proloquo2go, predictable etc. (large number of apps)	Direct Touch, switch access dependent on the communication app	no	Depends on the app and iPad case
Windows Tablet: surface pro /Laptop	Windows 10	Tablet 11 inches	Communicator 5 Grid 3 Compass Snap + Core First	Touch (tablet)/ mouse/joystick, 1 and 2 switch scanning (require switch interface or eye mobile Plus), eye gaze (requires PC Eye mini hardware or eye mobile plus bracket)	Yes with addition of EyeR dongle or Eye Mobile Plus bracket which has IR capability	no

SET BC Supported AAC Software/Apps Comparison Chart

Communication Software and Apps	Operating System hardware options	Page sets available and size	Symbol Sets available	Visual scenes	Editing Software	Environmental Controls	Keyboard options/word prediction	Access Options
Nuvoice (prentke romich) v. 2.10	Windows 10/ Can only be used on a Prentke Romich "Accent" AAC device.	-Page based: 4, 8, 15, 45 loc. -Core Scanner: 1, 8, 23, 25, 108, 190 location -Lamp WFL: 84 location -Unity: 4, 8, 15, 28, 36, 45, 60, 84, 144 location -(Picture) Word Power PCS/PRC: 45, 60, 84 location -Essence	-Unity (minspeak) -PCS -Symbol stix (add-on) -Clarity (add-on) - photos	Yes, with buttons and hotspots	Yes/ free windows Nuvoice Pass https://www.prentrom.com/support/PASS/download-nuvoice-pass-software	Yes. Pre made page sets or can be created and customized	-Qwerty, -ABC, AEIOU, - scanning, -60 Qwerty, -frequency, -keypad keyboards all with word prediction option	Touch, mouse, 1 and 2 switch scanning , auditory scanning, eye gaze
Communicator 5 (Tobii Dynavox) 5.4.12 (Pseudo operating system)	Windows 10 -surface pro tablet -Tobii I- series	-Emerging: sono primo 32 loc. w/ visual scenes -Symbol based: My words, My visual scenes, sonoflex 35 loc, sonoLexis 32, 54, 108 loc. LiterAACy 35, 54, 108 loc. Word power (add-on) -Text: sono scribe (letters, words phrases), sono key (social media) -Activity based/envir. control	-PCS, symbolstix -photos	Yes, with buttons and hotspots	Yes/trial or evaluation version, free full version for SLPs https://www.tobii-dynavox.com/en-US/products/mytobiidynavox/tobii-dynavox-for-professionals/?redirect=true	Yes. Pre made page sets or can be created and customized	-Several customizable Qwerty, abc, scan, pointer, dwell free. -Word and phrase prediction.	Touch, mouse/alt. joystick, eye gaze, 1 and 2 switch scanning
Grid 3 (3.0.38) (Pseudo Operating system)	Windows 10 -surface pro tablet -Tobii I series	Supercore 12, 24, 54 location Wordpower 60 , 100 Vocabulary for Life 42 loc. Symbol talker A, B, C, D, 24, 35, 48, 63 loc. B-Line: 99 location (motor planning) PODD 15, 60 location	Photos, Widget literacy symbols, PCS, symbol stix, (dependent on grid set)	Yes, button (hotspots ?)	Trial version	Yes, premade page sets	-Qwerty and ABC keyboards -word prediction available	Touch, eyegaze, mouse, 1 and 2 switch scanning

Communication Software and Apps	Operating System options	Page sets available and size	Symbol Sets available	Visual scenes	Editing Software	Environmental Controls	Keyboard options/word prediction	Access Options
Snap + Core First (Tobii Dynavox) v. 1.6.4487.0	Windows 10 -surface pro tablet -Tobii I series iOS/ipad (not currently available)	Communication page Grid size 1 to 80 location	Photos and PCS	no	Trial version or free version for SLPs	Yes, pre-made pages and can be edited with limited actions.	Qwerty, ABC and AZERTY (European French) -no word prediction for smaller grid sizes	Touch, eyegaze, switch scanning (limited)
Compass/ Compass w/ PODD (Tobii Dynavox) iApp v. 2.5.0.13481	Windows 10/Tablet windows/Tobii I- series iOS/iPad,	5 levels of ability, core first, word power, PODD, NavBar, Masterpage, text-based, word power, Gateway	-PCS -Photos	No but has scenes with pop up grids	Yes/trial version and free full version for SLPs	Some in windows devices	-Qwerty, word prediction	Touch, mouse, eye gaze, 1 and 2 switch, (dependent on hardware)
Lamp Words for Life iApp v. 2.18.0	-Windows 10/ Accent - iOS/iPad, -Android/Nova chat	84 location one hit 84 location transition 84 location full	-Unity -PCS available on Accent	no	Only for windows Nuvoice PASS (Accent)	Only for Windows version on Accent	-Qwerty keyboard -keyguards available for each device.	Dependent on individual hardware systems
NovaChat app v.	Android/Nova Chat device	Basic 4 loc., Multichat 15 Student, adolescent, adult, My Core, Quick Chat 4, 8, 12 child, 8, 12 adult, spelling, Vocab PC 12 Adult, adolescent, Word (chat) power 20 simply, 24, 30, 42 basic, 42, 48, 60 basic, 60, 80, 108	Symbolstix, PCS, Photos	yes	Yes/chat editor v. 2.12.0	no	-customizable QWERTY -Keyguards available	-Touch, -1 and 2 switch, -auditory scanning
GoTalk Now iApp v. 4.11.23	iOS/iPad	1, 4, 9, 16, 25,36 (customizable)	Photos, comes w/ Imagine, In-app purchase of symbolstix Widget & PCS	yes	no	no	customizable	Touch, auditory cues, in-app switch control (1 and 2 switches)

Communication Software and Apps	Operating System options	Page sets available and size	Symbol Sets available	Visual scenes	Editing Software	Environmental Controls	Keyboard options/word prediction	Access Options
TouchChat iApp v. 2.18.0	iOS/iPad	Basic 4 loc.; Multichat 15 (Student, adolescent, adult); My Core; Quick Chat (4, 8, 12 child, 8, 12 adult); spelling; Vocab PC (12 Adult, adolescent); Communication Journey: Aphasia; Word power (20 simply, 25 touch and scan, 42 basic, 42, 48, 60 basic, 60, 80, 108, 108 w/ keyboard); i Essence; Talk about AAC; Pixon 60; Gateway (child 12, 20, 40, 66, teen 20, 30, 40, 66 teen/adult, 66 text, Pro) French: (Free) 4 page sets http://touchchatapp.com/page-sets	-Symbolstix, photos, unity/minspeak symbols with Pixon page sets	yes	Yes/chat editor v. 2.12.0 http://touchchatapp.com/support/software-updates	no	-Customizable QWERTY, -Keyguards available depending on case	-Touch, -switch control through iOS (not in-App)
Proloquo2Go iApp v. 6.1	iOS/iPad	Core: 9 – 64 location Basic: 9 – 36 location	Symbolstix, photos	no	no	no	Qwerty, A - Z	In-App switch control (1 or 2 switch), touch, auditory cues
Sono Flex iApp v.2.0.0	iOS/iPad (Available as full version in Communicator 5)	32 location page set	Symbolstix, photos	no	no	no	ABC keyboard	Touch or iOS switch control
Scene & Heard iApp v. 3.12	iOS/iPad	-prebuilt scenes as examples for different purposes.	-photos, 12,000 widget symbols, drawing option, photo library	yes	no	no	n/a	Touch, in-app switch access: 1 or 2 switch