

BIGmack / LITTLEmack - Switches

General Description

The BIGmack is a flat 5-inch single message device that allows a user to record a message of up to 2 minutes in length.

The LITTLEmack is an angled 2.5-inch single message device that records a message of up to 75 seconds in length.

Supports students who:

are nonverbal and may have cognitive, physical, or visual challenges. The BIGmack is recommended for individuals who require a larger target.

Vendor website

AbeNet



Where Can I Learn How to Use It?

- [Getting Started with the Bigmack Speech Device](#)
AbleNet's website for, getting started guides, videos and other product support links for the BIGmack.
- [Getting Started with the LITTLEmack Speech Device \(AbleNet\)](#)
AbleNet's support hub for product questions, getting started guides, instructions for use, and other downloads on the LittleMack
- [Big Mack / Little Mack](#)
SET-BC training video on how to use the Bigmack and LITTLEmack.

Where Can I Get Ideas on How to Use It with Students?

- [101+ ideas for using the BIGmack – CENMAC](#)
pdf on Spectronix site.
- [30 ideas for you and your BIGmack! from Liberator by Prentke Romich](#)
powerpoint slides

What If I Want to Know Even More?

- [106 Ideas for using a BIGmack from TalkSense](#)
A comprehensive list of suggestions for using the BIGmack and LITTLEmack including best practices.
- [Using a LittleMack or BigMack With a Toy from Adaptive Tech Solutions](#)
This video shows how to use a BIGmack or LITTLEmack to activate a toy. Includes important info that the toy will activate only with a recorded message.